



## General Search Results--Full Record

Article 4 of 4

☒ UNMARK

[PREVIOUS](#)

[SUMMARY](#)

[HOLDINGS](#)

[FIND RELATED RECORDS](#)

[Explanation](#)

### MICRO-BRIDGES FOR SITTING-DROP CRYSTALLIZATIONS HARLOS K JOURNAL OF APPLIED CRYSTALLOGRAPHY 25: 536-538 Part 4 AUG 1 1992

Document type: Note   Language: English   Cited References: 7   Times Cited: 35

#### Abstract:

Micro-bridges, small devices in the shape of a bridge, have been designed to carry out **sitting-drop crystallizations** on a microscale. Micro-bridges have an indentation in the top of the bridge to contain the protein **drop** during a **crystallization** experiment and to prevent the droplet from spreading over a large area. The special feature of micro-bridges is that they fit into the **wells** of standard tissue-culture multiwell **plates**.

#### KeyWords Plus:

PROTEIN CRYSTALLIZATION

#### Addresses:

HARLOS K, OXFORD CTR MOLEC SCI, MOLEC BIOPHYS LAB, REX RICHARDS BLDG, S PARKS RD, OXFORD OX1 3QU, ENGLAND

#### Publisher:

MUNKSGAARD INT PUBL LTD, COPENHAGEN

#### IDS Number:

JH630

#### ISSN:

0021-8898

Article 4 of 4

[PREVIOUS](#)

[SUMMARY](#)

[Acceptable Use Policy](#)

Copyright © 2002 [Institute for Scientific Information](#)

HOME

Products

Stuff

Customer Service

Technical Support

About Us

Products -> Crystallization Plates and Accessories -> Micro-Bridges™

Fast Find... ▾

## Micro-Bridges™



Micro-Bridges™, small devices in the shape of a bridge, are designed to carry out sitting drop vapor diffusion crystallization when placed in a VDX or Linbro Plate. A single Micro-Bridge fits neatly into the reservoir of a standard 24 well VDX Plate for a sitting drop crystallization experiment. Once placed inside the wells, Micro-Bridges are stable and there is no need to stick them to the wells with grease or adhesive. It is therefore possible to transfer them to other wells during or after a crystallization experiment. Why sitting drop? Placing the droplet in the indentation greatly reduces the risk of losing the protein by accident. Crystallization can be carried out in the presence of detergents and organic solvents which are compatible with polystyrene (such as MPD, iso-propanol and ethanol). The protein drop is less affected by condensation problems. Soaking and seeding experiments can be carried out easily. Crystals can be transported more securely. Larger drop volumes can be used. Micro-Bridges have a round indentation in the top surface of the bridge which holds the sample droplet during a crystallization experiment and prevents the droplet from spreading over a large area. Made from polystyrene, these parts are highly transparent and suitable for most crystallizations. The surface of the indentation is highly polished to facilitate the visual inspection of the drops under a microscope. The maximum drop volume for the Micro-Bridge is 35 microliters. Reservoirs can be sealed with plain cover slides and vacuum grease or clear sealing tape.

### ⌚ Ordering Info

⌚ Catalog Number: HR3-310  
Description: Micro-Bridges  
Quantity: 100 pack sampler

Catalog Number: HR3-312  
Description: Micro-Bridges  
Quantity: 400 pack

[ Support | Customer Service ]

Hampton Research  
27632 El Lazo Road  
Suite 100  
Laguna Niguel, CA 92677-3913

Tel: 949-425-1321  
Fax: 949-425-1611

Copyright 2001  
All Rights Reserved.

HOME

Products

Stuff

Customer Service

Technical Support

About Us

Products -> Crystallization Plates and Accessories -> Micro-Bridges Polypropylene™

Fast Find...

## Micro-Bridges Polypropylene™



Same great Micro-Bridges, new material. Manufactured from clarified polypropylene, these Micro-Bridges are resistant to most organic solvents, and are especially useful with crystallization experiments that involve detergents and other hydrophobic reagents. The polypropylene Micro-Bridges resist drop spreading typically observed when using reagents such as detergents, volatile organics such as iso-propanol and ethanol, and non-volatile organics such as MPD. These Micro-Bridges also resist acetone, dioxane, acetonitrile, 2,2,2 trifluoroethanol, and other aggressive organic solvents.

### Ordering Info

**Catalog Number: HR3-340**  
**Description: Micro-Bridges Polypropylene**  
**Quantity: 100 pack sampler**

**Catalog Number: HR3-342**  
**Description: Micro-Bridges Polypropylene**  
**Quantity: 400 pack**

[ Support | Customer Service ]

Hampton Research  
27632 El Lazo Road  
Suite 100  
Laguna Niguel, CA 92677-3913

Tel: 949-425-1321  
Fax: 949-425-1611

Copyright 2001  
All Rights Reserved.